visomat

medic pro

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Table of contents

A	Intended use	17
B	Safety instructions	18
C	Operating the unit	
	1. Unit description	19
	2. Important instructions for use	21
	3. Measuring blood pressure	22
D	Technical information	
	1. Customer service	24
	2. Technical data	24
	3. Original spare parts and accessories	25
	4. Metrological inspection	26
	5. Explanation of symbols	26
	6. Disposal	26
E	Maintaining the unit	27
F	Warranty	28

Thank you for choosing the visomat medic pro upper-arm blood pressure monitor (also referred to in the following as the unit).

Intended use

A

This device is a mechanically operating blood pressure monitor for the auscultatoric measurement of the systolic and diastolic pressure in the upper arm. The measured blood pressure values can be read on the large scale of the manometer.

There are no restrictions as to the age, weight, physical or psychological condition of the patients when using the blood pressure measurement.

It is ideally suited for professional use in hospitals, in private practices, and in emergency services, providing safe monitoring, control and diagnosis of the blood pressure values.

Safety instructions

- Before putting on the cuff, please check if the arm circumference is in the range specified on the cuff band (see page 19). More cuff sizes are shown on page 25 accessories.
- The unit must only be operated with original parts. The warranty will be invalidated if the unit is damaged by non-approved accessories!
- Never open or modify the unit or the cuff - this is a medical device and may only be opened by authorised specialists. If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution.
- Make sure that the cuff and the air hose are not bent.
- Please comply with the specified operating conditions for measurement. See page 24 Technical data.
- The unit must be used in accordance with the procedures contained in these instructions and must not be used for other purposes.
- Shockproof in accordance with DIN EN ISO 81060-1:2012.
- Before every application, check if the pointer is still resting at zero point ± 3 mmHg. The pointer must move freely across the entire scale. If the device has been damaged, it needs to be technically inspected. See page 26 metrological inspection.


Operating the unit

1. Unit description


Pressure gauge with display dial

Cuff

Markers on the blood pressure monitor

- | | |
|--|---|
| 1 Product name | 5 Measuring unit for blood pressure (millimetre mercury column) |
| 2 TÜV Süd Product Service, Munich | 6 Country of origin |
| 3 Measurement range 0-300 mmHg | 7 Manufacturer |
| 4 Shockproof in accordance with DIN EN ISO 81060-1:2012. | 8 Serial number |
| | 9 Scale: 1 scale line = 2 mmHg |

2. Important instructions for use

- To keep the device ready for use, you have to connect the tube of the cuff with the corresponding connection of the manometer.
 - Blood pressure measurements of the patient should be performed under relaxed and calm conditions.
 - The size of the cuff needs to be selected to fit to the upper arm.
 - Make sure that there are no clothes obstructing the blood flow of the upper arm (e.g., by pulling up the sleeves).
 - Position the arm and palm facing upwards.
 - In order to rule out differences between left and right and to obtain comparable measurement results, it is important always to take blood pressure measurements on the same arm.
 - Please pay attention to the following recommendations of the EN ISO 81060-1:
 - In pregnant patients, K5* should be applied. On rare occasions – if Korotkoff sounds are still audible after completely venting the cuff – it is advisable to use K4** to determine the diastolic blood pressure. In any case, the correct measurement method should be determined in agreement with a doctor.
 - To measure blood pressure levels of children between the age of 3 and 12 years, K4** of the Korotkoff sounds should be used to determine the diastole.
- *K5 is the phase in which the Korotkoff sounds cannot be heard any more.
 **K4 is the point at which the frequencies heard through the stethoscope change from a clear beating to a muffled thumping sound.

3. Measuring blood pressure

- Expose the upper arm.
- Pull the cuff over the upper arm so that the lower edge is positioned 2-3 cm above the elbow joint, the tube should be leading downwards.
- Place the marker "arteria" above the artery of the upper arm. The control line for assessing the arm circumference must remain within the marker for the arm circumference. If this is not the case, choose a cuff fitting to the arm circumference of the patient.
- Close the air extraction vault by tightening the air drain screw.
- If blood-pressure values are not known, inflate the cuff to approx. 200 mmHg. However, if you can clearly hear thumping sounds in the stethoscope, continue inflating until no more sounds can be heard.
- Slowly unscrew the air drain screw, thereby allowing the air to escape from the cuff.
- The pointer of the manometer will slowly run back to its neutral position. By turning the air drain screw, the speed at which the air is vented can be adjusted. (Recommendation of the WHO: 2-3 mmHg/second) Now observe exactly the needle excursion on the dial and listen out for the noises which you can hear through the stethoscope.
- On hearing the first clearly thumping sound, read off the systolic blood-pressure value on the dial and make a note of the reading (Fig. 1).
- Listen out for rhythmic knocking in your stethoscope. When the sounds stop or at the moment of the last audible sound, read off


Fig. 1

the diastolic blood-pressure value on the dial and make a note of this reading too (Fig. 2).

- After determining the diastole, open the air drain screw to release the remaining air from the cuff. Remove the cuff after the air is completely drained.


Fig. 2

1. Customer service

The unit may only be repaired by the manufacturer or by an expressly authorised organisation. Please contact:

UEBE Medical GmbH
Bgm.-Kuhn-Str. 22
97900 Kilsheim, Germany
info@uebe.com
www.uebe.com

Please only send the device together with the cuff.

2. Technical data

Model:	visomat medic pro
Unit type:	Aneroid blood pressure monitor
Weight:	Approx. 200 g without cuff
Measurement procedure:	Auscultatory measurement
Measurement and pressure-display range	0-300 mmHg
Measurement precision:	Cuff pressure: ± 3 mmHg
Serial number:	The unit bears a serial number which provides clear identification. First two digits: year of manufacture, third digit: quarter, remaining digits: consecutive numbering.
Cuff:	Available from an arm circumference of 8-55 cm, see page 25 Accessories
Operating conditions:	Ambient temperature 10 to 40 °C Relative humidity up to 85%, non-condensing
Storage and transportation conditions:	Ambient temperature -20 to +60 °C Relative humidity 10% to 95%, non-condensing

Recommendations for storage:	dry and at room temperature
Time to adjust from minimum/maximum storage temperature to operating temperature:	2 hours
Deflation:	Manual

3. Original spare parts and accessories

The following original spare parts and accessories are available from specialist dealers:

- Cuff type 1 (for arm circumference 8-15 cm)
Part no. 2700007, PZN-09786255
- Cuff type 2 (for arm circumference 14-21,5 cm)
Part no. 2700005, PZN-09786261
- Cuff type 3 (for arm circumference 18-25 cm)
Part no. 2700003, PZN-09786284
- Cuff type 4 (for arm circumference 26-33 cm)
Part no. 2700001, PZN-09786290
- Cuff type 5 (for arm circumference 34-43 cm)
Part no. 2700002, PZN-09786309
- Cuff type XXL (for arm circumference 43-55 cm)
Part no. 2700008, PZN-09786315

For further spare parts, please contact Customer Service (see page 24).

Subject to technical modifications.

4. Metrological inspection

The metrological properties need to be tested at least every 2 years when using the device in a commercial setting. This can be implemented either by UEBE Medical GmbH, an authority responsible for metrology, or authorised maintenance services. For this, please observe your national provisions.

Upon request, responsible authorities or authorised maintenance services receive a "Test instruction for metrological inspection" from the manufacturer.

Important: No modifications, e.g. opening the unit, may be made to this unit without the manufacturer's prior permission.

Please only submit the device for metrological inspection together with the cuff and the instructions for use.

5. Explanation of symbols

This product complies with the Council Directive 93/42/EC from 5 September 2007 regarding medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Manufacturer

Reference number = item number

Observe the safety instructions

6. Disposal

At the end of its service life, you may dispose of the device with normal household waste.

- The unit contains sensitive parts and must be protected against strong variations in temperature, air humidity, dust and direct sunlight.
- Make sure when storing the unit that no heavy objects are left lying on the unit or on the cuff and that the air hose is not kinked.
- The unit is not water-proof.
- If the unit becomes dirty during use, use a soft, dry cloth to clean the unit. Do not use benzene, thinners or other strong solvents.
- To disinfect the device and the diaphragm of the stethoscope, we advise using Mikrozid Sensitive Liquid to avoid damage to the synthetic material and the sight glass.
- Do not scrub or machine-wash the cuff. If the cuff gets dirty during use, use a synthetic cleaner and softly rub the surface.
- After removing the air pocket, the case of the cuff can be washed at 30 °C. To reinsert the air pocket, feed the air tube through the air tube opening from the inside. Make sure that the air pocket with the longer end is pulled through the cuff first. Please ensure that the air pocket is placed into the cuff without any folds and creases.
- To disinfect the cuff, we recommend using a disinfection spray with 70% isopropyl.
- If the unit has been opened, it must be subjected to a metrological inspection conducted by an authorised organisation.
- Make sure that liquid cannot get into the air hose. Dry thoroughly.
- If the device has been stored at temperatures below 10°C or above 40°C, allow it to stand at room temperature for at least 2 hours before using it.

The device has been manufactured and tested with great care. However, in the unlikely event of a defect being detected after delivery, we provide warranty in accordance with the following terms and conditions:

During the warranty period of 5 years from the date of purchase we reserve the right either to repair any such defect at our expense or to supply a perfect replacement unit.

The warranty does not cover normal wear and tear of wearing parts or damage caused by failure to observe the instructions for use, improper handling (e.g. breakage) and/or disassembly of the device by the buyer. Furthermore, no claims for damages against us are substantiated by the warranty.

Warranty claims can only be advanced in the warranty period and by presenting proof of purchase. In the event of a warranty claim, the unit must be sent to the following address together with the proof of purchase and a description of the complaint:

UEBE Medical GmbH
Service-Center
Bgm.-Kuhn-Str. 20
97900 Kulsheim
Germany

The cost of returning the unit to our factory shall be borne by the sender. UEBE shall refuse to accept return deliveries that have not been paid for by the sender.

The statutory claims and rights of the buyer against the seller (claims for defect, manufacturer's liability, for example) are not restricted by this warranty.

Please note: In the event of a warranty claim it is essential to attach the proof of purchase.

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© Copyright 2020 UEBE Medical GmbH