# visomat

## comfort form

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## Table of contents

<b>A</b>	<b>Intended use</b>	40
<b>B</b>	<b>Safety instructions</b>	41
<b>C</b>	<b>Operating the device</b>	
	1. Device description	44
	2. Display	45
	3. Important usage information	46
	4. Initial operation of the device	46
	5. Inserting/replacing batteries, battery symbol	47
	6. Setting the date, time and sounds	48
	7. Activating and deactivating the button lock	50
	8. Attaching the cuff	50
	9. Body position during the measurement	51
	10. Selecting the user	52
	11. Measuring blood pressure	52
	12. Motion control	56
	13. Classification of the measured values (WHO)	56
	14. Irregular pulse	57
	15. Using the memory	57
<b>D</b>	<b>What you should know about blood pressure</b>	
	1. Systolic and diastolic blood pressure value	60
	2. Why you measure different values	60
	3. Why should you measure blood pressure regularly?	60
<b>E</b>	<b>Technical information</b>	
	1. Error messages	61
	2. Customer service	62
	3. Technical data	62
	4. Genuine spare parts and accessories	63
	5. Metrological inspection	64
	6. Explanation of symbols	64
	7. Electromagnetic compatibility (EMC)	65

## Table of contents

<b>F</b>	<b>Cleaning and care</b>	69
<b>G</b>	<b>Warranty</b>	70

Thank you for choosing the visomat comfort form upper arm blood pressure monitor (hereinafter also referred to as the "device").

These instructions for use are intended to help the user operate the device safely and efficiently. They must be kept with the product and passed on as necessary. The unit must be used in accordance with the procedures contained in these instructions for use and must not be used for any other purposes.

It is important that you read all the instructions carefully before using the device. Please pay particular attention to the chapter "Important usage information" on page 46.

The intended purpose of this device is to enable persons over 12 years of age to automatically measure their own systolic blood pressure, diastolic blood pressure and pulse at their upper arm. The device can also be used by health care professionals.

### Contraindications

The device is not suitable for patients undergoing dialysis treatment or taking anticoagulants, platelet aggregation inhibitors or steroids. Using the cuff in these circumstances may cause internal bleeding.

- Wash your hands after every measurement. If the device is used by more than one person, wash your hands before and after every use.
- Taking blood pressure measurements on children requires specialist knowledge. Consult your doctor if you want to measure a child's blood pressure. Never use this device on a newborn baby or an infant.
- The device is not suitable for use by pregnant women.
- Do not under any circumstances place the cuff over a wound or an aneurysm or on an arm with an arteriovenous shunt. Risk of injury!
- When the cuff is inflated, any supply via an intravascular access point (infusion) or other medical monitoring devices could be interrupted.
- Do not take a blood pressure measurement on the same side as a mastectomy with simultaneous removal of the axillary lymph nodes.
- The device does not have category AP/APG classification. It is not suitable for use in the presence of flammable anaesthetics containing air, oxygen or nitrogen oxide.
- The device is not suitable for use in close proximity to electrosurgical instruments or magnetic resonance imaging equipment.
- The device cannot be used as a replacement for professional ECG monitoring equipment for heart rate monitoring.
- The device contains small parts that children could swallow. The air hose poses a strangulation hazard. Do not leave the device unattended with children.
- Store the device out of the reach of children and pets.
- Ensure that the air hose is not kinked. A kinked air hose may prevent the cuff from deflating, thus interrupting the blood flow in the arm for too long.


## Safety instructions

- Please wait a few minutes between successive measurements, otherwise the blood flow in the arm is interrupted for too long and injuries could occur.
- If the measurement indicates an irregular pulse (arrhythmia) (see page 57), consult with a doctor before evaluating the results.
- The readings measured using this device do not provide a diagnosis. They are no substitute for visiting a doctor. Do not rely solely on the result obtained; always take the patient's other symptoms into account too. In case of doubt, please call a doctor or paramedics to help.
- Do not assess the measurement results yourself. Only health care professionals are qualified to interpret blood pressure measurements. Discuss any unusual measured values with your doctor. Never change the dose of medicines prescribed by your doctor!
- The unit contains sensitive parts and must be protected from extreme temperature fluctuations, humidity, shocks, dust and direct sunlight.
- Do not knock or drop the device. Avoid exposing the device to strong vibrations.
- The unit must only be operated with the cuff and accessories provided (see "Genuine spare parts and accessories" on page 63). Using different cuffs or accessories can lead to incorrect measurement results. The warranty will be invalidated if the unit is damaged by third-party accessories!
- When using the device with a mains adapter, where possible, use the visomat mains adapter, which has been specially tested for medical devices (see „Genuine spare parts and accessories“ on page 63). External adapters connected to medical electrical equipment via the DC input must comply with the relevant IEC or ISO standards. In addition, all configurations must meet the requirements for medical electrical systems (see IEC 60601-1-1 or Section 16 of the 3rd edition of IEC 60601-1 or 60601-1-2). Anyone connecting external adapters to medical electrical equipment is configuring a medical system and is therefore responsible for ensuring that the system meets the require-

## Safety instructions


- ments for medical electrical systems. Please note that local laws take precedence over the above requirements. If in doubt, please contact visomat customer service (see page 62).
- Using this device in the vicinity of mobile phones, microwaves or other devices with strong electromagnetic fields can cause malfunctions and inaccurate measurements.
- Never open or modify the device or the cuff (except for replacing the batteries). If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution
- The intended user is the patient. The device must be used in accordance with the information on measurement, battery replacement and maintenance contained in these instructions for use. It must not be used for any purpose other than those described in these instructions for use.
- To avoid inaccurate readings, please adhere to the intended operating and storage conditions. See "Technical data" on page 62.
- The inflating and measuring procedure can be interrupted at any time by pressing the Start/Stop button or by removing the cuff. The device then stops the inflation procedure and deflates the cuff.
- Before use, check that the device is operating reliably and is in proper condition. Never use the device if it is damaged. Doing so can lead to injuries and permanent damage or inaccurate measurements.
- When using the power adapter, always make sure that the adapter and cable are not damaged.
- Never attempt to service or repair the device yourself. If the device malfunctions or is faulty, please refer to the troubleshooting guide starting on page 61, or contact customer service (see page 62).
- Any serious adverse effects linked to the product must be reported to the manufacturer and the relevant authorities in the member state in which the user/patient resides.

## 1. Device description

- 1** LCD display
- 2** Colour scale for WHO classification
- 3** Connecting socket for mains adapter
- 4** Start/Stop button
- 5** Battery compartment
- 6** Memory button
- 7** Set button (SET)
- 8** User selection button
- 9** Cuff connection
- 10** Air hose
- 11** Artery marking
- 12** Measurement arrow for arm circumference
- 13** Marking for arm circumference

## 2. Display

- 1** User 1
- 2** Date/time
- 3** Battery display
- 4** Motion control
- 5** SYS = systolic (top blood pressure number)
- 6** Blood pressure graph
- 7** WHO classification (page 56)
- 8** DIA = diastolic (bottom blood pressure number)
- 9** Irregular pulse
- 10** Pulse signal display
- 11** Memory number
- 12** PUL = calculated pulse frequency per minute
- 13** Average value memory symbol
- 14** Triple measurement symbol
- 15** Cuff position check
- 16** User 2

### 3. Important usage information

- Do not consume alcohol, caffeine or smoke for at least one hour before measuring.
- Rest for at least 10 minutes before taking the measurement. Depending on the intensity of any preceding exertion, you may even need to rest for up to an hour.
- Expose the upper arm. Clothing must not hinder the blood flow to or from the arm, as this will affect the blood pressure at the measuring point and may make the readings inaccurate.
- Do not move or talk during measuring.
- Breath slowly and deeply. Do not hold your breath.
- In order to obtain comparable measurement values, always take the measurements under comparable conditions where possible; for example, always at the same time and in the same place.
- In order to exclude side differences and to obtain comparable measurement results, it is important to always take the measurement on the same arm. Ask your doctor which side it is advisable for you to measure.
- Blood pressure is not a fixed value. It can go up or down by 20–40 mmHg in patients within a few minutes.

### 4. Initial operation of the device

Insert the batteries provided into the device.

If the device is to be operated using mains power, the cable connector of the mains adapter (not supplied) must be plugged into the socket on the right-hand side of the device. The batteries are automatically switched off.

Where possible, please use the visomat mains adapter, which has been specially tested for medical devices. See “Genuine spare parts and accessories” on page 63. When using a different mains adapter, always observe the technical specifications for the mains adapter (see “Technical data” on page 62).

### 5. Inserting/replacing batteries, battery symbol

- Remove the battery compartment cover on the underside of the device.
- Remove the old batteries from the device and insert the new ones. Please ensure correct polarity (markings in battery compartment).
- Close the battery compartment again. Make sure that the battery cover is securely closed.

Only use long-life alkaline batteries (4 x 1.5 V, size AA/LR6).

When replacing batteries, always replace all batteries simultaneously. Do not mix new and old batteries or batteries of different types. Do not use batteries beyond the specified expiry date.

After the “Battery low” symbol has appeared on the display for the first time during inflation (Figure 1), you can carry out approx. 30 more measurements, but you will need to have replacement batteries ready.

Figure 1

If the “Battery flat” symbol flashes in the display (Figure 2), no further measurements can be taken. Switch off the device and replace the batteries.


Figure 2

Once the batteries have been replaced, the readings will still be stored in the memory, but the date and time will need to be reset. To do this, follow the steps below.

If the device will not be operated for more than three months, remove the batteries. Leaking battery fluid can damage the device. If leaked battery fluid comes into contact with eyes, rinse the affected eye immediately with plenty of water. Seek immediate medical attention.

## 6. Setting the date, time and sounds

- Once the batteries have been inserted during initial commissioning of the device or after replacing the batteries, the date and time must be set again. The year flashes on the display (Figure 1).


Figure 1

- Press the User button to reduce the value by 1 or the Memory button to increase the value by 1. To quickly increase or reduce the value, hold down the relevant button. Set the correct year and confirm the entry by pressing the SET button.
- You are now able to set the month (Figure 2). Press the User button to reduce the value by 1 or the Memory button to increase the value by 1. To quickly increase or reduce the value, hold down the relevant button. Set the correct month and confirm the entry by pressing the SET button.


Figure 2

- Repeat these steps to set the day (Figure 2), the hours and the minutes (Figure 3).

Figure 3

- After setting the date and time, you are taken to the sound settings. Here you can switch the beeps that can be heard during the measurement on and off.

Figure 4

- Press the User button or the Memory button to switch between SP on (speaker on, Figure 4) and SP off (speaker off, Figure 5).

- Confirm your selection by pressing the SET button.

Figure 5

### Changing the date, time, and sound settings

- Press and hold the SET button until the year flashes on the display.
- Set the values as described on page 48.

### Displaying the date and time

You can also display the date and time when the device is switched off by pressing the SET button. The date and time are shown alternately on the display. The display will switch off again after approx. 60 seconds.

## 7. Activating and deactivating the button lock

The device is equipped with a button lock to prevent it from being switched on unintentionally, e.g. when you are travelling and carrying the device in your luggage.

To activate the button lock, press and hold the User button and the Memory button simultaneously for 3 seconds when the device is switched off. "Loc on" appears on the display (Figure 1).

Figure 1

Now if you press a button on the device, the display shows "Loc on" (Figure 1). However, the device does not switch on.

Figure 2

To deactivate the button lock, press and hold the User button and the Memory button simultaneously for 3 seconds. "Loc on" appears on the display (Figure 1). Keep the buttons held down until "Loc off" appears on the display (Figure 2). The button lock is now released and you can use the device again as normal.

## 8. Attaching the cuff

Before applying the cuff, check that the circumference of your upper arm is within the range specified on the cuff. Wrap a tape measure around the centre of your upper arm when it is relaxed to measure its circumference. An incorrect cuff size could result in inaccurate readings.

- Expose your upper arm.
- Place the cuff around your left upper arm. The lower edge of the cuff should be 2-3 cm above the elbow (Figure 1). The artery marking on the cuff (ARTERY) must be


Figure 1

placed over the artery that runs along the inside of the arm.

- Now pull the free end of the cuff tight and close the Velcro fastener (Figure 2). The cuff should not be too tight. You should be able to insert two fingers between the arm and the cuff.


Figure 2

- Check that the measurement arrow for the arm circumference on the closed cuff is within the "Marking for arm circumference" on the cuff edge (Figure 3).


Figure 3

- Connect the plug of the cuff to the cuff connection on the left-hand side of the device. Ensure that the plug is fully inserted into the device. Do not force the plug into the hole.

## 9. Body position during the measurement

- Sit comfortably at a table (preferably at the height of a dining table, not a coffee table!). Sit with your back against the backrest of the chair
- Place the arm with the cuff on the table so that it is loosely stretched out. Hold it very still during measurement and do not talk. The palm of your hand should be facing upwards.
- Put your feet on the floor and do not cross your legs.

## 10. Selecting the user

The device is designed for two users. The values for the two users are stored separately and can also be evaluated separately.

You may also use the two user memories to evaluate measurements taken in the morning and those taken in the evening separately. To do this, carry out the measurement in the morning using User 1 and in the evening using User 2.

Before taking the measurement, check that the correct user is shown on the display.

- To change the user, press the User button. The currently set user appears on the display. Press the User button again to switch between User 1 and User 2.
- Now carry out the measurement.

## 11. Measuring blood pressure

You can choose to perform either one single measurement or one triple measurement, in which three automatically controlled single measurements are performed successively with brief intervals between them.

To perform a single measurement, press the Start/Stop button briefly. To start a triple measurement, press and hold the Start/Stop button for at least three seconds.

You can cancel the measurement at any time by pressing the Start/Stop button.

### Performing a single measurement

- Switch on the device by pressing the Start/Stop button. The full display is shown (Figure 1).


Figure 1

- After the device has completed its calibration against the ambient air pressure, the automatically controlled inflation process starts. The cuff is inflated to the pressure required for measurement.

- Shortly after the inflation procedure begins, a check is carried out to ensure that the cuff has not been applied too loosely. If the cuff is applied correctly, the  symbol is displayed (Figure 2).

Figure 2

If the symbol is not displayed by the time the pulse signal starts, the cuff is too loose. In this case, you should repeat the measurement after you have applied the cuff correctly (see page 50).

- The actual measurement process then begins as the cuff deflates (Figure 3). Once the device has detected your pulse, the  pulse rate icon flashes.

Figure 3

- While the measurement is being performed, you have the opportunity to change the user again. To do this, press the User button to switch between User 1 and User 2. Once the measurement is complete, it is no longer possible to change the user, as the measured values have already been saved in the respective user memory by this point.
- Once the measurement has been taken, the cuff automatically deflates. The determined systolic, diastolic and pulse values are shown on the display (Figure 4).


Figure 4

## Operating the device

- A bar chart appears on the right of the display showing the measurement results for the last 7 days. Here you can track your blood pressure measurement history (Figure 5). The right-hand bar shows the current measurement. The bars for the previous 7 days appear to the left. If several measurements were performed on the same day, the average value of those measurements is displayed. If no measurements were taken on a particular day, the column remains empty.
- The device switches off automatically after approx. 60 seconds. You can also switch the device off before this by pressing the Start/Stop button.


Figure 5

### Performing a triple measurement

- Press and hold the Start/Stop button for at least three seconds to switch on the device and to start the triple measurement. The full display is shown (Figure 6).
- The first measurement will then begin. The triple measurement symbol (Ⓢ) is shown on the display with the number of the measurement (1, 2, 3) to the right of it (Figure 7).
- The continued sequence of measurement is described under "Performing a single measurement" on page 52.


Figure 6

Figure 7

## Operating the device

- At the end of the measurement, the result is shown in the display. A bar chart appears on the right of the display showing the results of the individual measurements. The current measurement is displayed in the right-hand column. (Figure 8).
- There is an interval of two minutes to prevent the blood flow in the arm from being interrupted for too long. A clock counts down from 120 seconds to 0 next to the triple measurement symbol (Ⓢ) (Figure 8).
- After the 120 seconds elapses, the second measurement starts automatically. After the second measurement is complete, there is another interval of 120 seconds before the third measurement.
- The current result (second measurement) is now displayed in the right-hand column of the bar chart. Column 1 shows the result of the first measurement (Figure 9).
- Once the third measurement is complete, the average values are calculated and shown on the display (Figure 10). The average result of the triple measurement is displayed in the right-hand column of the bar chart. Columns 1–3 show the results of the three individual measurements.


Figure 8


Figure 9

Figure 10

The average values of the triple measurement are calculated from the values of the last two measurements. The values of the first measurement are disregarded.


The device switches off automatically 10 seconds after the end of the measurement. You can also switch the device off before this by pressing the Start/Stop button.

## 12. Motion control

Movements during the measurement may result in incorrect measurement values. The motion control registers stronger movements and these are indicated by the ▼ symbol in the display. Other disruptive factors such as speaking, coughing or cuff noises can cause the symbol to be displayed. If this symbol appears, please repeat the measurement and avoid the disruptive factors.

## 13. Classification of the measured values (WHO)

The device classifies the measured blood pressure values according to the guidelines of the World Health Organization (WHO) and the International Society of Hypertension (ISH) of 1999. The triangle on the right of the display and the colour scale printed on the display edge show you the classification after each measurement:

Classification	Systolic pressure = top value mmHg	Diastolic pressure = bottom value mmHg
Severe hypertension (grade 3)	≥ 180	≥ 110
Moderate hypertension (grade 2)	160–179	100–109
Mild hypertension (grade 1)	140–159	90–99
High-normal	130–139	85–89
Normal	120–129	80–84
Optimal	< 120	< 80

Depending on age, weight and general health, the blood pressure values can be different. Only a doctor can determine the right blood pressure range for you and assess whether your blood pressure has reached a dangerous level for you. Discuss your blood pressure values with your doctor. **Never change the dose of medicines prescribed by your doctor!**

## 14. Irregular pulse

The ▼ symbol indicates that definite pulse irregularities have been detected during the measurement. In this case, the measurement result may deviate from your normal blood pressure. Repeat the measurement.

If the symbol appears, this is usually no cause for concern. However, if the symbol appears on a regular basis (e.g. several times in a week when blood pressure is measured daily) or more often than usual, consult a doctor for advice.

The device must not be used in place of a cardiac exam; it is only intended for early detection of an irregular pulse.

## 15. Using the memory

The measured results are automatically saved in the memory. The device has two memories (one per user) for storing measured values. Each memory can store up to 120 readings and the average value.

### Data retrieval

- Select the user whose data you want to retrieve (see page 52).
- With the device switched off, press the Memory button.

## Operating the device

- If the most recent measurement performed was a single measurement, the average value from the last three single measurements is calculated and displayed with "AVG" (Figure 1).


Figure 1

- A bar chart appears on the right of the display showing the average values for the last 7 weeks. Here you can track your blood pressure measurement history over a longer period of time. You can use the SET button to switch between the display of the last 7 weeks (Figure 2) and the display of the last 7 days (Figure 3).


Figure 2

Figure 3

- Press the Memory button again to display the most recently measured values (Figure 4). The saved value display also alternates between showing the date and time of the measurement.


Figure 4

- If the last measurement performed was a triple measurement, the average value of the triple measurement is indicated with "AVG" and the triple measurement symbol. The average value of a triple measurement is calculated from the last two


Figure 5

measurements. The result of the first measurement is disregarded during the calculation. The bar chart shows the average value of the triple measurement in the right-hand column. Columns 1–3 show the results of the individual measurements (Figure 5).

Use the User button and the Memory button to scroll up and down through the memory.

## Operating the device

The device switches off automatically approx. 30 seconds after a button was last pressed.

### Deleting data

- Select the user whose data you want to delete (see page 52).
- With the device switched off, press the Memory button to access the memory.
- To delete individual measured values, recall the measured value you wish to delete.
- Now press and hold the SET button for at least 4 seconds until the display flashes before being replaced by dashes (Figure 6).
- To delete the entire memory, retrieve the respective average value (AVG.). Now press and hold the SET button for at least 4 seconds until the average value flashes before being replaced by dashes (Figure 6).


Figure 6

## D What you should know about blood pressure

### 1. Systolic and diastolic blood pressure value

The cardiovascular system has the important task of supplying all organs and tissues in the body with sufficient blood, and of removing metabolites. The heart does this by contracting and expanding at a regular rhythm of around 60-80 times per minute. The pressure created by the blood flow on the artery walls when your heart contracts is referred to as systole. The pressure in the subsequent relaxation phase when the heart fills with blood again is referred to as diastole. In your daily measurement, you will measure both values.

### 2. Why you measure different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. Affected by mental, physical and environmental influences, it varies all the time and never remains constant.

Reasons for fluctuating blood pressure values can be: moving, speaking, eating, consuming alcohol or nicotine, nervousness, inner tension, emotions, room temperature, recent urination or bowel movements, environmental influences such as movements and sounds, etc. Even changes in the weather and climate can have an effect on your blood pressure.

This explains why values measured at the doctor are often higher than those you obtain at home in your usual environment.

### 3. Why should you measure blood pressure regularly?

Even the time of day has an effect on your blood pressure. The values are at their highest during the day. In the course of the afternoon and in the evening, they drop slightly. They are low while you are sleeping, but rise again relatively quickly once you get up.

One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is only possible when measurements are taken regularly. Discuss the measured values with your doctor.

## Technical information

### 1. Error messages

Error	Possible cause	Remedy
Er 1 shown on the display	Sensor error or faulty pressure pump	If the fault occurs several times, contact the customer service centre.
Er-1 shown on the display	Pulse measurement under 40 or over 199.	No device error. Discuss any values outside of the measuring range with your doctor.
Er 2 shown on the display	Device could not detect a pulse or could not calculate blood pressure values.	Check the hose line and air plug for blockages. Repeat the measurement. If the error occurs again, please contact the customer service centre.
Er 3 shown on the display	Abnormal measured values	Repeat the measurement. If the fault occurs several times, contact the customer service centre.
Er 4 shown on the display	The cuff takes a long time to inflate. The cuff is not sitting correctly or the connection to the tube is not airtight.	Fasten the cuff correctly and check that the air plug has been fully inserted into the device. Repeat the measurement.
Er 5 shown on the display	The air hose is kinked or the cuff is fastened too tightly.	Unkink the air hose. Correctly fasten the cuff. Repeat the measurement.
Er 6 shown on the display	The device has detected extensive pressure fluctuations.	Repeat the measurement. Do not move. Do not speak.
Er 7 shown on the display	The pressure measured by the sensor is over the limit value.	Repeat the measurement. If the fault occurs several times, contact the customer service centre.
Er 8 shown on the display	The limit is incorrect or the device does not have a limit.	Repeat the measurement. If the fault occurs several times, contact the customer service centre.
Err is shown on the display, device switches off	Arm movement as the cuff was being pumped	Repeat the measurement. Do not move.
	Talking during the measurement	Repeat the measurement. Do not speak.
No display after the device is switched on	Batteries inserted incorrectly	Check position of batteries
	Batteries flat	Replace batteries

Error	Possible cause	Remedy
The cuff is not inflated	Cuff not connected correctly	Check that the air connector is positioned correctly and repeat the measurement
	Air connector or cuff damaged	Replace cuff (see "Genuine spare parts and accessories" on page 63)

## 2. Customer service

Device repairs may only be performed by the manufacturer or an expressly authorised body. Please contact:

UEBE Medical GmbH  
 Bgm.-Kuhn-Str. 22  
 97900 Kulsheim,  
 Germany  
 info@uebe.com, www.uebe.com

## 3. Technical data

Model:	visomat comfort form
Size:	L = 138 mm x W = 100 mm x H = 45.7 mm
Weight:	303 g without batteries
Display:	LCD display (liquid crystal display) 130 x 70 mm
Memory:	2 memories each with 120 measured values (automatic storage) and average value (AVG)
Measurement method:	Oscillometric measurement
Pressure sensor:	resistive
Pressure display range:	0–295 mmHg
Measuring range:	Diastolic: 40–220 mmHg Systolic: 60–260 mmHg Pulse: 40–199 pulses/min.
Measurement accuracy:	Blood pressure measurement (cuff pressure): $\pm 3$ mmHg, pulse rate: $\pm 5\%$

Measurement resolution:	1 mmHg
Power supply:	4 x 1.5 V AA LR6 alkaline manganese batteries
	Optional: Mains power supply, class II, input: 100–240 V AC 50/60 Hz, output 5 V DC, 1 A. Compliance with IEC 60950-1, IEC 60601-1
Operating conditions:	Ambient temperature 5 to 40°C, relative humidity 15 to 90%, air pressure 700 to 1060 hPa
Conditions for storage and transport:	Ambient temperature -20 to 60°C, relative humidity 10 to 93%, air pressure 700 to 1060 hPa
IP rating:	IP 21: Protected against solid foreign bodies with a diameter of 12.5 mm and above, protection against dripping water.
Protection against electric shock:	Internal power supply, applied part type BF (cuff)
Expected service life:	5 years
Operating mode:	Continuous operation

## 4. Genuine spare parts and accessories

The following genuine spare parts or accessories are available through specialist retailers:

Special shell-shaped cuff type PF1 for arm circumference 22–42 cm  
 Art. no. 2403801, PZN-19925410

D-Ring cuff type VWR3 for arm circumference 22–42 cm  
 Art. no. 2402801, PZN-19925427

D-Ring cuff type WXL2 for arm circumference 32–52 cm  
 Art. no. 2407801, PZN-19925404

visomat mains adapter USB-C/Micro-USB  
 Art. no. 2400030, PZN-19779500

Subject to technical amendments.

## 5. Metrological inspection

Generally, a metrological inspection is recommended every 2 years. However, professional users in Germany are obligated to do so in accordance with the "Medical Products Operator Ordinance" (Medizinprodukte-Betreiberverordnung).

The inspection can be performed either by UEBE Medical GmbH, an authority responsible for metrology or an authorised maintenance service. Please refer to your national regulations.

Upon request, responsible authorities or authorised maintenance services will receive a "test instruction for metrological inspection" from the manufacturer.

Please only submit the device for metrological inspection together with the cuff and the instructions for use.

Caution: No modifications such as opening or repairing the device may be made without the manufacturer's permission (with the exception of replacing the batteries).

## 6. Explanation of symbols

**CE** 0123 This product complies with Regulation (EU) 2017/745 of the European Parliament and the Council dated 5 April 2017 regarding medical products. It bears the marking CE 0123 (TÜV SÜD Product Service GmbH).

**MD** Medical product

**REF** Reference number = item number

**SN** Device serial number

**UDI** Unique product identifier

 Consult the instructions for use

 Caution

 Manufacturer

 Temperature limits

 Humidity limits

 Store in a dry place

 Degree of protection against electric shock: TYPE BF

 Direct current

 Protection class II (double insulation)

**IP21** Protected against solid foreign bodies with a diameter of 12.5 mm and above, protection against dripping water.

 Technical equipment and batteries do not belong in household waste. They must be disposed of at appropriate collection and disposal points.

## 7. Electromagnetic compatibility (EMC)

The device is medical electrical equipment and is subject to special EMC precautions, which must be published in the operating instructions.

Portable and mobile RF communications equipment may affect the equipment. Using the device with non-approved accessories can adversely affect the device and change the EMC.

The device should not be used immediately adjacent to or between other electrical devices.

The device meets the EMC requirements set out in international standard IEC60601-1-2. The requirements are met under the conditions described below.

#### Guidance and manufacturer's declaration – electromagnetic emissions

Emission measurements	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic currents according to IEC 61000-3-2	Class A
Voltage fluctuations/flicker according to IEC 61000-3-3	Met

#### Guidance and manufacturer's declaration – electromagnetic immunity

Test frequency	Modulation	Immunity test level
30 kHz	Continuous wave	8 A/m
134.2 kHz	Pulse modulation <sup>b</sup> 2.1 kHz	65 A/m <sup>c</sup>
13.56 kHz	Pulse modulation <sup>b</sup> 50 kHz	75 A/m <sup>c</sup>

a This test is only applicable to medical electrical devices and systems intended for use in home health care.

b The carrier must be modulated with a square-wave signal with a duty factor of 50%.  
c RMS value, before the modulation is applied

#### Guidance and manufacturer's declaration – electromagnetic immunity

Immunity tests	IEC 60601-1-2 test level	Compliance level
Electrostatic discharge (ESD) according to IEC 61000-4-2	±8 kV contact discharge ±2 kV, ±4 kV, ±8 kV, ±15 kV air discharge	±8 kV contact discharge ±2 kV, ±4 kV, ±8 kV, ±15 kV air discharge
Electrical fast transient/bursts according to IEC 61000-4-4	±2 kV for power lines ±1 kV for input and output lines 100 kHz repetition rate	±2 kV for power lines
Surges according to IEC 61000-4-5	±0.5 kV voltage phase conductor–phase conductor ±2 kV voltage phase conductor–earth ±1 kV voltage phase conductor–phase conductor	±0.5 kV voltage phase conductor–phase conductor ±1 kV voltage phase conductor–phase conductor
Voltage dips, short interruptions, and voltage variations according to IEC 61000-4-11	0% for 0.5 periods at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% for 1 period and 70% for 25/30 periods. Single phase at 0° 0% for 300 periods	0% for 0.5 periods at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% for 1 period and 70% for 25/30 periods. Single phase at 0° 0% for 300 periods

Immunity tests	IEC 60601-1-2 test level	Compliance level
Power frequency magnetic fields according to IEC 61000-4-8	30 A/m 50 Hz/60 Hz	30 A/m 50 Hz/60 Hz
RF interference conducted through cables according to IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms (in ISM and amateur radio frequency bands) 80% AM at 1 kHz	3 Vrms 150 kHz to 80 MHz 6 Vrms (in ISM and amateur radio frequency bands) 80% AM at 1 kHz
Radiated RF interference according to IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz
Radiated fields in close proximity according to IEC 61000-4-39	30 kHz: 8 A/m 134.2 kHz: 65 A/m 13.56 MHz: 7.5 A/m	30 kHz: 8 A/m 134.2 kHz: 65 A/m 13.56 MHz: 7.5 A/m

NOTE: U<sub>i</sub> is the AC mains voltage before the test level is applied

#### Guidance and manufacturer's declaration – electromagnetic immunity

Radiated RF interference according to IEC 61000-4-3 – test specifications for the immunity of enclosures to high-frequency wireless communication devices				
Test frequency MHz	Frequency band <sup>a</sup> MHz	Radio service <sup>a</sup>	Modulation <sup>a</sup>	Immunity test level V/m.
385	380 to 390	TETRA 400	Pulse modulation <sup>b</sup> 18 Hz	27
450	430 to 470	GMRS 460, FRS 460	FM <sup>c</sup> ± 5 kHz hub 1 kHz sinus	28
710	704 to 787	LTE band 13, 17	Pulse modulation <sup>b</sup> 217 Hz	9
745				
780				
810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE band 5	Pulse modulation <sup>b</sup> 18 Hz	28
870				
930				
1720	1700 to 1990	GSM 1800; CDMA 1900; GSM 1900; DECT;	Pulse modulation <sup>b</sup> 217 Hz	28
1845				
1970				
2450	2400 to 2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450 LTE band 7	Pulse modulation <sup>b</sup> 217 Hz	28

Test frequency MHz	Frequency band <sup>a</sup> MHz	Radio service <sup>a</sup>	Modulation <sup>a</sup>	Immunity test level V/m.
5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation <sup>b</sup> 217 Hz	9
5500				
5785				
NOTE: If necessary, the distance between the transmitting antenna and the medical electrical device or system can be reduced to 1 m in order to achieve the immunity test level. The 1 m test distance is permitted according to IEC 61000-4-3.				
<p>a For some radio services, only the frequencies for the radio connection from the mobile communication device to the base station (uplink) have been included in the table.</p> <p>b The carrier must be modulated with a square-wave signal with a duty factor of 50%.</p> <p>c As an alternative to frequency modulation (FM), pulse modulation with a 50% duty cycle at 18 Hz can be used, as this would represent the worst case, if not the actual modulation.</p>				

- Keep the device clean. Check for cleanliness after use. Please use a soft, dry cloth for cleaning. Do not use benzene, thinners or other strong solvents.
- The cuff contains a sensitive current of compressed air. Handle the cuff carefully and avoid any type of strain due to twisting or kinking. Keep the cuff away from sharp or pointed objects.
- The cuff may absorb sweat and other fluids. Inspect the cuff for marks or discolourations after each use. To clean the device, carefully wipe the surface using a damp soapy cloth. 70% isopropyl alcohol can also be used if necessary. Do not wash the device in the dishwasher or washing machine. Do not immerse the device in water.
- Take care when storing the device. Ensure that no heavy objects are resting on the device or the cuff and that the air hose is not kinked. Do not wind the air hose too tightly.
- To separate the cuff from the main device, please do not pull the air hose, but grip the air connector and gently pull it off.

The device has been manufactured and tested with all due care. Nevertheless, in the event of defects upon delivery, we provide a warranty under the following conditions:

During the warranty period of 5 years from the date of purchase, we remedy such defects at our discretion and at our expense in our plant through repair work or replacement delivery of a defect-free device.

The warranty does not cover normal wear and tear of wearing parts or damage caused by failure to observe the instructions for use, improper handling (e.g. unsuitable power sources, breakage, leaking batteries) and/or disassembly of the device by the buyer. Furthermore, the warranty does not constitute grounds for asserting claims for damages against us.

Warranty claims can only be asserted during the warranty period and upon presentation of the proof of purchase. In the case of a warranty claim, the device along with the proof of purchase and a description of the complaint must be sent to:

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The cost of returning the device to our factory shall be borne by the sender. Complaints that are sent back without prepayment will not be accepted.

The statutory claims and rights of the buyer against the seller (for example, defect claims, producer liability) are not restricted by this warranty.

Please note: In case of warranty, please attach the proof of purchase.

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© Copyright 2025 UEBE Medical GmbH